



FIELD CHANGE ORDER

Service

Issued by : DMC Hamburg

Released : G. Kramm **signed**

Publication No. : 4512 980 56521

Reference No. : 00 135 006

Date : February 2000

Product Group : 742

*This FCO replaces the issue
(98.0) of September 1998*



OPTIMUS 50 / 65 / 80 RAD / RF

APPLIES TO:

All OPTIMUS 50 / 65 / 80 RAD / RF generators

9890 000 02001

9890 000 00871

9890 000 00661

9890 000 00851/2

9890 000 00881

9890 000 00862

all Optimus generators up to and including

serial numbers

970218 - 970220 970223 - 970340

970361 - 970442 970444 - 970446

970447 - 970763

971002 - 971130

972601 - 972632

980007 - 980445

982001 - 982186

982601 - 982700 982998 + 982999

TITLE:

Upgrade of Release 3.x Optimus generators to Release 3.3

LIST OF PAGES & DRAWINGS:

1 ... 3 (00.0)

INTRODUCTION:

- Symptom** :
- 00B5 desk lookup
 - 00CQ during system turn on
 - no APR backups possible via APRMAN or XRG-Scope
 - release of fuse F3 on mA-control PCB
 - error combinations
 - 00CD - 08CF
 - 00CX - 08CY
 - 00CF - 08CY
 - erros
 - 00CD, 00CY, 00CQ, 08CY
 - errors during the switch - on phase
 - 00XI - 00M3
- Cause** : Interferences and software instabilities can lead to different problems.
- Remedy** : Release upgrade to release 3.3 including additional interference suppression activities.

MANPOWER / TIME TO COMPLETE:

1 engineer 3.5 hours

TOOLS & TEST EQUIPMENT:

- 3 + 4 mm Allan keys
- 10 mm open-end wrench
- low watt soldering iron with **fine tip**
- service PC with hardkey and mouse
- flash light
- PLCC extraction tool AMP 822154-1
2422 487 89772
- standard tools

MODIFICATION KIT / PARTS REQUIRED:

For all generators mentioned under APPLIES TO, the following parts can be ordered at SSA Hamburg.

4512 104 39893 Upgrade kit OPTIMUS rel. 3.3

comprising:

4512 113 20524	Eprom Human Interface
4512 104 39901	EZX45W_X52 CAN Termination plug
2022 333 00024	MKP-KO 33N/10% 1000V
4512 104 90531	Cabel extension SYSTEM CAN EZ X43
0322 042 00015	Wire for N5:23 - N5:26
4512 116 02203	"Service Disk A"
4512 984 08482	SERV.MAN.UPGR.KIT

Note:

The modification kit will be available Free Of Charge till the end March 1999 at the spare part centre in Hamburg. Please refer to the generator serial number when ordering the modification kit.

PROCEDURE:

- Exchange of firmware Eprom in Human Interface.
- Software Upload Release 3.3 into CU PCB
- Installation of a short circuit in Rotor Control HS unit.
- Installation of CAN extension cable.

Detailed instruction manual supplied with the upgrade kit.

PARTS DISPOSAL:

All parts are to be disposed of in a safe way in accordance with local Safety Regulations.

DOCUMENTATION:

File this FCO and the Upgrade instruction manual in binder: OPTIMUS under Tab 8

